



FACE ROLLER

HOT & MICROCURRENT MASSAGER



FACE ROLLER

ULTIMA EMS FACE ROLLER

Advanced Home Aftercare for Ultima Face Treatments

- Bridges the gap between in-clinic treatments and daily home care
- Designed specifically as an aftercare extension of Ultima Face
- Supports visible results between sessions
- Easy to introduce, easy to sell, easy to use



ULTIMA

WHY AFTERCARE MATTERS

Professional results need daily support

- In-clinic treatments initiate change
- Daily at-home care maintains and enhances results
- Structured aftercare increases client satisfaction
- Better aftercare = stronger retention

ULTIMA



WHY THE FACE ROLLER IS ESSENTIAL FOR ULTIMA AFTERCARE

- Helps maintain firmness and facial definition
- Supports lymphatic drainage and reduces puffiness
- Keeps facial tissues active between treatments
- Strengthens the professional-to-home treatment journey
- Encourages long-term engagement with Ultima



KEY SELLING POINTS

- Combines EMS microcurrent, warm therapy and cold therapy
- Gentle, non-invasive and suitable for daily use
- No learning curve – intuitive for all clients
- Compact, elegant and travel-friendly
- Designed to complement Ultima Face protocols



TECHNOLOGY: EMS MICROCURRENT

Gentle muscle activation for facial support

- Low-frequency EMS (8.3 Hz)
- Supports facial muscle tone and contour
- Helps maintain lifting effects between sessions
- Comfortable and safe for regular use

Microcurrent Stimulates skin

Anti-wrinkle and skin reguvenation,Maintains youthful-looking skin



TECHNOLOGY: WARM THERAPY

Comfortable warmth for circulation and relaxation


- Supports microcirculation
- Helps relax facial muscles
- Prepares skin for better product absorption
- Ideal before or during skincare application

40±3°C Hot Compress

More Than Tighten Skin, Import Nutrients Effectively

40±3°C temperature is slightly warmer than human body, which gently and effectively opens pores.

EMS micro-current promote muscle movement, which greatly improves absorption of skin care products.



TECHNOLOGY: COLD THERAPY

Cooling sensation for soothing and de-puffing

- Helps reduce puffiness and fluid retention
- Supports lymphatic flow
- Refreshes tired-looking skin
- Perfect for morning or post-treatment use

Roll Out Tighter And Glossy Skin

Persist in using it with skin care products for more than 40 days, restore firming and glossy skin



WHY IT'S BETTER THAN BASIC FACE ROLLERS

- Active technologies vs. passive rolling only
- Warm & cold functionality in one device
- Designed specifically for facial aftercare
- Premium experience aligned with clinical environments



PRICE & MARGINS

Recommended Retail Price (RRP): 160 EUR

- Attractive distributor margins
- Competitive premium positioning
- Excellent impulse upsell after treatments



TECHNICAL INFORMATION

- EMS frequency: 8.3 Hz
- Battery: 3.7 V / 1100 mAh
- Charging time: approx. 4 hours
- Water resistance: IPX4
- Weight: 96 g
- Packaging: Premium gift box





THANK YOU!

WWW.ULTIMASOLUTION.COM